

New ! Stand Alone Tray & Load Control Labeling Systems



Item Description		Hospital Name
DBL PUNCTURE LAP TUBAL TRAY		
PROCESS STEAM	STORE IN 2nd. Floor	TECH T.LEE
		1095A

COLORS AVAILABLE - IN STOCK !

Load Control Electronic Systems



HOSPITAL NAME / DEPARTMENT		
STERILIZER #	LOAD #	TECH I.D.
04	07	JC1
STERILE  V-PRO® UNLESS OPENED, DAMAGED or WET		
DATE	TIME	COUNTER
12/25/18	10:55AM	01



800.679.7836
WWW.S1MEDICAL.COM

AUTHORIZED DISTRIBUTORS:



ITEM IDENTIFICATION AT ITS BEST TRAY AND BUNDLE LABELING SYSTEM

Create, save & print on demand custom labels
with these complete Stand Alone Systems
No PC Required.

COMPACT + MOBILE CONVENIENT + EFFICIENT

- Perfect Quality Label Print
- No Need to Stock Printed Labels
- Automatic Date & Time Print
- Pre or Post Sterilization Labeling



SM-4 and SM-1 Stand Alone Printers Communication & Tray Labeling

Synthetic Labels for all item / tray identification
Pre & Post Sterilization Labeling



SM-8 Stand Alone Printers Load Control Labeling

Paper or Synthetic Labels
Pre and Post Sterilization Labeling

WHICH SYSTEM IS BEST FOR YOU

SM SYSTEM MODELS

				
	Z70-SM1KT	Z70-SM4KT	Z70-SM8KT	Z70-SM8KD
ITEM STORAGE SPACE	1,000 With Mem. Upgrade	4,500 STANDARD	1,000 STANDARD	1,000 STANDARD
TRAY IDENTIFICATION (SYNTHETIC LABELS)	✓	✓	✗	✗
AUTOMATIC DATE AND TIME PRINT	✓ OPTIONAL	✓ OPTIONAL	✓	✓
STEAM, GAS & STERRAD VALIDATED LABELS	✓	✓	✓	✗
COMPATIBLE WITH ALL TRAY TRACKING SYSTEMS	✓	✓	✓	✗
PRINT EXTRA LARGE FONTS	✓	✓	✓	✓
SMALL FOOT-PRINT - COUNTER SPACE SAVER	✓	✗	✓	✓
PRE-STERILIZATION LABEL CAPABLE	✓	✓	✓	✗
FULL ITEM DESCRIPTION/SERVICE/ TECH, ETC.	✓	✓	✓	✓
AUTOMATIC TRAY, BUNDLE or ITEM COUNTER	✗	✗	✓	✓
BARCODE SCANNER READY (OPTIONAL)	✓	✓	✓	✓
FULL KEYBOARD w/ 4 LINE DISPLAY READY	✓	✓	✓	✓
PRE-PROGRAMMED & READY-TO-USE OUT OF BOX	✓	✓	✓	✓

SM8 SYSTEMS - LOAD CONTROL LABELING

Labels available in paper and synthetic for labeling before or after sterilization

HOSPITAL AND DEPARTMENT IDENTIFICATION	HOSPITAL NAME / DEPARTMENT			TECHNICIAN IDENTIFICATION
STERILIZER AND LOAD DATA	STERILIZER # 04	LOAD # 08	TECH I.D.: JR02	BARCODE INFORMATION
COLOR PRINT EVENT RELATED STATUS & DISCLAIMER	STERILE V-PRO UNLESS DAMAGED, OPENED or WET			STERILIZATION VALIDATED SYSTEMS MATERIAL
AUTOMATIC DATE & TIME	12/25/18	TIME 8:45 AM	COUNTER 01	TRAY, BUNDLE OR ITEM COUNTER



SM8t Printer




SM8d Printer

AVAILABLE PRINTS



CATALOG NO.: **L51-A55EY**
SM8t-SYNTHETIC



CATALOG NO.: **L51-A55EB**
SM8t-SYNTHETIC



CATALOG NO.: **L51-P50DW**
SM8d-PAPER DIRECT

TRAY IDENTIFICATION - SAMPLE LABEL FORMATS

COLORS AVAILABLE MULTIPLE SIZES

HOSPITAL NAME / DEPARTMENT		DATE: 12/25/18 TIME: 10:45 AM
MIS Modular Cutting Blocks Knee		
PROCESS V-PRO	SERVICE CARDIO	TECH G.M.
1-6001		

Hospital Name	4758-Z1
HAYHURST VASCULAR CLAMPS	
PROCESS STEAM	SERVICE VASCULAR

Hospital Name	787D
STRYKER HIP RESURFACING TRAY - 01	
PROCESS STERRAD	STORE IN ORTHOPEDIC

COMMUNICATION - SAMPLE LABEL LAYOUTS

COLORS AVAILABLE MULTIPLE SIZES

STERRAD
System Sterilized

STERILE
UNLESS DAMAGED, OPEN or WET
* CHECK CAREFULLY BEFORE USE *

LOANER EQUIPMENT

INDEFINITE SHELF LIFE
STERILE
UNLESS DAMAGED OPEN OR WET
CHECK CAREFULLY


NOTICE IMPLANTABLE ITEM

This Tray is **INCOMPLETE** Date: _____
Missing: _____


Select Size and Design Your Perfect Label !

Design over the label size that fits best your needs and fax or e-mail it to us at 813-972-0578 or sales@system1us.com. We will e-mail back to you a professional label design for you to see your label in perfect detail. Note which fields are to be large and bold for best results.



STOCK LABEL SIZES colors: Yellow, Blue, Pink, Green




SIZE: 4" x 2" CAT.# **L51-A170** (COLORS)




SIZE: 4" x 1.5" CAT.# **L51-A205** (COLORS)



SIZE: 2" x 1.25"
CAT.# **L51-A210** (COLORS)



SIZE: 3" x 5" CAT.# **L51-A110** (WHITE)



SIZE: 3" x 2"
CAT.# **L51-A32** (COLORS)

CATALOG NUMBER

SM STAND ALONE SYSTEMS & ACCESORIES

Z70-SM4K

Printer ZT230T, keyboard, 20 pre-loaded label layouts, software & cables

Z70-SM1Kt

Printer GX420T, keyboard, 20 pre-loaded label layouts, software & cables and labels

Z70-SM8Kt

Printer TLP2824, keyboard, 20 pre-loaded label layouts, software & cables and labels

Z70-SM8Kd

Printer LP2824, keyboard, 20 pre-loaded label layouts, software, cables and labels